

BARDIA

BARDIA Silicone Elastomer Foley Catheter 3-Way

AKL 20805502873


SPESIFIKASI

- Diproduksi dengan metode double blending sehingga integritas kateter terjaga dan laju aliran urin lebih konsisten.
- Volume balon 30cc. Penggunaan maksimal 7 hari

STERILE EO

UKURAN	PCS/BOX	KODE
18 FR	10	BARD-18(3W)
20 FR	10	BARD-20(3W)

UKURAN	PCS/BOX	KODE
22 FR	10	BARD-22(3W)
24 FR	10	BARD-24(3W)



Sterile unless or damaged.
Warning: Do not use lubricants having oil base as they will damage the catheter.
Do not aspirate drainage funnel.
Single use only. For urologic use only.
Valve Type: Unidirectional. Do not use needles.